



MiniMed™ Mio™ 30













Infusion set • Cathéter • Infusionsset • Equipo de infusión • Infusieset • Set di infusione • Slangesett
Infusiosetti • Infusionsset • Infusionssæt •
Conjunto de infusão • Infuzijski set • Infuzní set •
Infúziós szerelék • Zestaw infuzyjny • Infúzyon seti •
Infúzna súprava • Σετ έγχυσης • Инфузионный
набор • 输注管路 • أداة التثريب • ט עירווי •




Front • Avant • Vorderseite • Cuerpo de frente • Voorkant • Davanti • Foran • Edessä • Fram • Forside • Frente • Spredaj • Čelní strana • elölnézet • Przód • Ön • Vpredu • Εμπρός • Живот • 正面 • الأمام • חיזית

Back • Arrière • Rückseite • Cuerpo de espalda • Achterkant • Dietro • Bak • Takana • Bak • Bagside • Costas • Zadaj • Zadní strana • hátulnézet • Tyl • Arka • Vzadu • Πίσω • Спина • 背面 • الخلف • גב




DISCONNECTING • DÉCONNEXION •
 TRENNEN • DESCONEXIÓN •
 LOSKOPPELEN • SCOLLEGAMENTO •
 KOBLE FRA • IRROTTAMINEN •
 BORTKOPPLING • FRAKOBLING •
 DESLIGAR • ODKLÓP • ODPOJENÍ •
 LEVÁLASZTÁS • ODŁĄCZANIE •
 BAĞLANTIYI KESME • ODPOJENIE •
 ΑΠΟΣΥΝΔΕΣΗ • ОТСОЕДИНЕНИЕ •
 断开连接 • الفصل • ניתוק

RECONNECTING • RECONNEXION •
 WIEDERANSCHLIESSEN •
 RECONEXIÓN • AANSLUITEN •
 RICOLLEGAMENTO • KOBLE TIL
 IGJEN • ΚΥΤΚΕΜΙΝΕΝ TAKAISIN •
 ÄTERANSLUTNING • GENTILKOBLING •
 VOLTAR A LIGAR • PONOVNÍ PRIKLOP •
 OPĚTOVNÉ PŘIPOJENÍ •
 ÚJRACSATLAKOZTATÁS • PONOWNE
 PODŁĄCZANIE • YENİDEN BAĞLAMA •
 OPĀTOVNÉ PRIPOJENIE •
 ΕΠΑΝΑΣΥΝΔΕΣΗ • ПОВТОРНОЕ
 ПОДСОЕДИНЕНИЕ • 重新连接 •
 חיבור מחדש • إعادة التوصيل

MiniMed™ Mio™ 30

EN See page 2

FR Consultez la page 6

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ES Consultar la página 14

NL Zie blz. 18

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NO Se side 26

FI Katso sivu 30

SV Se sidan 34

DA Se side 38

PT Consultar a página 42

SL Glejte stran 46

CS Viz strana 50

HU Lásd: 54. oldal

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TR Bkz. sayfa 62

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Indications for use

The MiniMed Mio 30 infusion set is indicated for the subcutaneous infusion of insulin from an infusion pump. The infusion set is indicated for single use.

Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

Description

MiniMed Mio 30 is an angled infusion set with an integrated insertion device. The built-in legs support a 30° insertion angle. MiniMed Mio 30 is delivered ready for use. Both MiniMed Mio 30 and the separate tubing unit are sterile.

Contraindications

The infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous (IV) infusion. Do not use the infusion set for infusion of blood or blood products.

Warnings

- Read all instructions for use before inserting the infusion set. Failure to follow instructions may result in pain or injury
- The infusion set is only sterile and non-pyrogenic if the package is undamaged and unopened. Do

not use the infusion set if the package is open or damaged. Ensure sterility by checking that the sterile paper and tamper-proof seal are not damaged.

- If using the infusion set for the first time, do so in the presence of a healthcare professional.
- Inaccurate medication delivery, site irritation, or infection may result from improper insertion and maintenance of the infusion site.
- Do not reinsert the introducer needle into the infusion set. Reinsertion may cause damage or tearing of the soft cannula, which may lead to unpredictable medication flow.
- Do not reuse the infusion set. Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.
- Carefully remove the needle guard before inserting the infusion set. Do not bend the needle prior to insertion. Do not use the infusion set if the needle is bent or has been damaged.
- If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.

Precautions

- Never try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.
- Do not leave air in the infusion set or the tubing. Make sure to fill the tubing completely. For filling instructions, consult the pump manufacturer's instructions for use.
- The infusion set is a single-use device which must be disposed of after use. Do not clean or resterilize.
- Check the infusion set frequently to ensure that the soft cannula remains firmly in place. Since the cannula is soft, it will not cause any pain if it slips out, and this can happen without notice. Replace with a new infusion set at a new site if the infusion set is not in place. The soft cannula must always be completely inserted to receive the full amount of medication.
- Do not put disinfectants, perfumes, alcohol, or deodorants on the infusion set, as these may impact the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose, or if the soft cannula is displaced from its original site. If the tape is curled, do not attempt to straighten it. Stretching or pulling on the tape may affect the cannula under the skin.

- Replace the infusion set every two to three days, or as directed by a healthcare professional.
 - If the infusion site becomes irritated or infected, replace the infusion set, and use a new site until the previous site has healed.
 - Use aseptic techniques when temporarily disconnecting the infusion set. Consult a healthcare professional for how to compensate for missed medication while the infusion set is disconnected.
 - To avoid accidental needle stick injury, never point a loaded insertion device toward any body part where insertion is not desired.
 - Check blood glucose one to three hours after insertion. Measure blood glucose regularly. Consult a healthcare professional for more information.
 - Do not change the infusion set just before bedtime, unless blood glucose can be checked one to three hours after insertion.
 - If blood glucose becomes inexplicably high, or an occlusion alarm occurs, check the infusion set for clogs or leaks. If in doubt, change the infusion set, because the soft cannula may be dislodged. Discuss a plan with a healthcare professional to lower blood glucose. Test blood glucose to make sure the problem is corrected.
 - Carefully monitor blood glucose levels when disconnected and after reconnecting the infusion set.
 - Disconnect the tubing with caution by pressing the grips. A hard pull on the tubing can result in damage to the infusion set. Ensure that the infusion set is properly in place when the tubing is fully disconnected.
- ### Serious incident
- If, during the use of this device, or as a result of its use a serious incident occurs, report the incident to the manufacturer and to the national authority.

Storage and disposal

- Store the infusion set in a cool, dry location at room temperature. Do not store the infusion set in direct sunlight or high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharps container according to local laws, to avoid the risk of needle stick injury.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

Warranty

For product warranty information, contact a local Medtronic support representative, or visit: www.medtronicdiabetes.com/warranty.

© 2020 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (TM) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.


MiniMed™ Mio™ 30

MiniMed Mio 30 instructions for use

Read all instructions carefully before using MiniMed Mio 30. Consult the pump user guide before connecting MiniMed Mio 30 to the pump. Follow proper hygiene procedures. If using MiniMed Mio 30 for the first time, do the first set-up in the presence of a healthcare professional.

Contents

- A. Lid
- B. Needle guard
- C. Soft cannula
- D. Adhesive tape
- E. Cannula housing
- F. Site connector
- G. Tubing connector
- H. Introducer needle
- I. Insertion device
- J. Sterile paper
- K. Finger grips – for holding MiniMed Mio 30

- L. Release buttons – press here when releasing MiniMed Mio 30
- M. Disconnect cover
- N. Protective cap for site connection

- O. Legs
- P. Spring with finger loop (inside)

Inserting MiniMed Mio 30

Unfold the front cover of this booklet to view images that are keyed to the steps below.

CAUTION: Check the tamper-proof seal and sterile paper. Do not use if the seal or paper has been broken.

1. Wash your hands before inserting the infusion set.
2. Select a recommended insertion site, shown in grey, and as indicated by a healthcare professional. Consult a healthcare professional for insertion site rotation.
3. Clean the area of insertion with a disinfectant as directed by a healthcare professional. Allow skin to air dry for at least 30 seconds before inserting the infusion set. If needed, remove hair around the insertion site to prevent improper attachment of the tape to the skin.
4. Remove the paper from the tubing. Connect the tubing to a filled reservoir.

IMPORTANT: Make sure there is no liquid inside the tubing connector or on top of the reservoir. If this occurs, start over with a new reservoir and infusion set.

5. Remove the protective cap of the site connector. Fill the tubing until insulin drips from the tip. Ensure there are no air bubbles in the tubing. Consult the pump manufacturer's instructions for use for proper filling procedures.

6. Pull the red tape to remove the seal.
7. Pull off the sterile paper. Do not touch the spring at this stage.
8. Hold the MiniMed Mio 30 with one hand and pull the lid off with the other.
CAUTION: Be careful not to bend or touch the MiniMed Mio 30 introducer needle. If the infusion set is not securely placed in the insertion device with the needle pointing straight ahead, pain or minor injury may occur during insertion.
9. Place your fingers on the finger grips on each side. With your other hand, use your finger to gently pull back the spring until you hear a “click”. The needle guard will lift up or peel off. Do not touch or bend the introducer needle. Check that the soft cannula does not exceed the introducer needle.
10. Be sure the adhesive tape is not stuck on the introducer needle.
11. **INSERTION:** Place your index finger on the top release button and your thumb on the bottom release button. Position MiniMed Mio 30 with its legs resting flat on the skin to ensure an insertion angle of 30°. Make sure to maintain a 30° angle while inserting and gently press the release buttons once to insert MiniMed Mio 30.
12. Make sure that the cannula stays in place. Gently use two fingers to hold down the tape while

removing the insertion device and introducer needle. Gently pull the insertion device straight back.

13. Secure the cannula by gently pressing down with one finger on the clear window of the tape. Carefully remove the protective liner under the cannula housing.
14. Remove the backing paper under the front end. Gently massage the adhesive tape onto the skin.
15. Put the lid back onto the MiniMed Mio 30 for disposal.
16. Gently place a finger on the cannula housing to hold it and insert the site connector into the infusion set using your other hand. Make sure the connector locks with a “click”. Fill the empty space in the cannula with 0.7 units (0.007 ml) of **U-100 insulin**.

Disconnecting

MiniMed Mio 30 allows you to temporarily disconnect your pump without having to change the infusion set.

- A. Gently place a finger in front of the cannula housing and squeeze both sides of the site connector. Pull the connector straight out.
- B. Insert the protective cap onto the site connector until you hear a “click”.
- C. Insert the disconnect cover into the cannula housing until you hear a “click”.

Reconnecting

- D. Remove the disconnect cover from the cannula housing.
- E. Remove the protective cap from the site connector. If necessary, fill the infusion set until insulin dispenses from the needle.
CAUTION: When filling, hold MiniMed Mio 30 with the needle pointing down.
- F. Gently place a finger in front of the cannula housing. Push the site connector into the cannula housing until you hear a “click”.

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Medtronic Ltd.
Medtronic Österreich GmbH
Unomedical Devices S.A. De C.V.